ARTHROSA^{MID}[®]

INSTRUCTIONS FOR USE

Arthrosamid®

en	4
cs	9
da	14
de	19
el	25
es	31
fi	36
fr	41
it	46
ko	51
no	56
pt	61
sv	66
tr	71

INSTRUCTIONS FOR USE

IMPORTANT

Please read all information carefully before using Arthrosamid®. Failure to properly follow instructions may result in improper functioning of the device which may lead to complications.

PRODUCT DESCRIPTION

Arthrosamid® is a non-absorbable, biocompatible, injectable, transparent hydrogel for intra-articular injection. Arthrosamid® consists of 2.5 % cross-linked polyacrylamide and 97.5 % non-pyrogenic water. Arthrosamid® is supplied in a pre-filled, sterile, 1 mL syringe, sealed with a Tip Cap in a Luer Lock fitting.

Arthrosamid® is intended to be injected intra-articularly with a sterile 21 G x 2 inches (0.8 x 50 mm) needle with Luer Lock fitting. A 3-part label with the LOT number, UDI and expiry date is available on the blister pack. Attach one label to the patient record to ensure product traceability.

Note: The scale on the syringe has no measuring function and is for guidance only.

HOW SUPPLIED

Arthrosamid® is sterilized by moist heat. Intended for single use only. Do not re-sterilize.

Upon removal from package, carefully inspect the product for signs of damage as this may compromise sterility and functionality. Do not use the product if the sterile barrier is damaged or unintentionally opened before use.

MODE OF ACTION

- Arthrosamid® integrates into the synovial tissue of the inner joint capsule and decreases joint stiffness
- Arthrosamid® forms a cushion-like membrane comprising the hydrogel incorporated in a matrix of fine tissue network covered by a synovial lining
- Arthrosamid® diminishes pain and improves function of the knee affected by osteoarthritis

The non-absorbable, non-bio-degradable and non-migratory characteristics of Arthrosamid® provide cushioning of the inner joint capsular tissue.

INDICATION

Arthrosamid® is intended to be used for symptomatic treatment of patients with knee osteoarthritis.

PATIENT TARGET GROUP

Arthrosamid® is intended for patients above the age of 18 years that have failed conservative treatments.

INTENDED USERS

Arthrosamid® must be administered by qualified physicians familiar with intra-articular injection procedures.

PERFORMANCE CHARACTERISTICS

Arthrosamid® is a permanent, non-degradable implant with an expected performance lifetime of 5 years.

CLINICAL BENEFITS TO BE EXPECTED

Arthrosamid® diminishes pain and stiffness and improves physical function of the knee and quality of life.

CONTRAINDICATIONS

Arthrosamid® should not be injected:

- If an active skin disease or infection is present at or near the injection site.

- If the joint is infected or severely inflamed.
- If the patient has previously received treatment with a different non-absorbable injectable/implant in the same knee.
- If the patient has previously received treatment with an absorbable intra-articular injectable (such as hyaluronic acid) that is not yet expected to be absorbed, according to manufacturer's information for the specific product.
- If the patient has received a knee alloplasty or has any foreign material in the knee.
- If the patient has undergone knee arthroscopy within the last 6 months.
- In patients with uncontrolled hemophilia or in patients with uncontrolled anti-coagulation treatment.

Do not subsequently inject a different non-absorbable implant in the injected knee.

WARNINGS

- Do not inject Arthrosamid® intravascularly, extra-articularly or in the tissue of the capsule.
- Do not inject corticosteroids along with Arthrosamid®.
- Do not inject any pharmaceuticals or biological substances along with Arthrosamid®.

PRECAUTIONS

Usual precautions associated with invasive joint procedures should be followed.

- Do not re-use Arthrosamid®. Re-use increases the risk of contamination and thereby increases the risk of infection.
- Do not re-sterilize Arthrosamid®.
- Do not use Arthrosamid® once expired.
- Do not inject Arthrosamid® into other parts of the body.
- Patients should avoid strenuous activity (e.g., tennis, jogging, or long walks) during the first weeks after the injection.

Arthrosamid® should be used with caution in patients:

- with acute or chronic infections in other parts of the body.
- with autoimmune disorders.
- in immunosuppressive therapy.
- with a compromised immune system.
- suffering from uncontrolled diabetes.
- undergoing major dental work or an invasive medical procedure.

Safety and effectiveness have not been established for:

- patients under 18 years of age.
- pregnant or lactating women.
- patients with BMI ≥ 35 kg/m².
- single doses higher than 6 mL.

As with any invasive joint procedure, there is a risk of infection when injecting Arthrosamid® and therefore prophylactic antibiotics must be administered prior to injection.

PATIENT INFORMATION

- The patient should be informed about indications, expected results, contraindications, warnings, precautions, and potential side effects.
- As with any invasive joint procedure it is recommended to avoid strenuous activity (e.g., tennis, jogging or long walks) during the first few weeks after the injection and begin rehabilitation activities, when possible.
- Advise patients that if they require a surgical or dental procedure in the future, they should tell the treating physician that they have a permanent implant and they should dis-

cuss the possible need for prophylactic antibiotic therapy with their treating physician.

- Patients should inform health care professionals of their Arthrosamid® implant for accurate future medical assessment.
- In case of complications, the patient should contact the injecting physician immediately for treatment.
- In the unlikely event that removal of Arthrosamid® is required, it can be removed surgically (synovectomy).
- Please provide the patient information leaflet and implant card to the patient.

SIDE EFFECTS

Common side effects:

Injection-related mild to moderate pain and/or mild swelling during the first weeks to months following the injection may occur. Formation or swelling of Baker's cysts may also occur.

Rare side effects:

Infections might occur at the injection site and must be treated immediately with antibiotics. In case of infection, the use of corticosteroids should be avoided as this may worsen the bacterial infection and complicate its treatment.

REPORTING OF INCIDENTS

All incidents must be reported to the local distributor or directly to Contura International A/S via Email to: complaints@contura.com.

If, during use of Arthrosamid®, you have reason to believe that a serious incident occurred, please report it to Contura International A/S and to your national authority.

DISPOSAL

Arthrosamid® may be disposed of as general waste. The needles used for the procedure must be disposed of as biohazardous materials (clinical waste) according to normal practice/local regulations.

STORAGE

The hydrogel must be kept away from sunlight. Do not freeze. The expiry date is indicated on the package.

METHOD OF ADMINISTRATION

PRE-PROCEDURE

Prophylactic antibiotics must be administered prior to injection with Arthrosamid®. The following combination of oral antibiotics administered 1-2 hours prior to Arthrosamid® injection is recommended:

For patients without penicillin allergy:

- Azithromycin 1 g plus Flucloxacillin 1 g.

For patients with known or suspected allergy towards penicillin:

- Azithromycin 1 g plus Clindamycin 600 mg or
- Azithromycin 1 g plus Moxifloxacin 400 mg.

These combinations of antibiotics achieve adequate tissue concentrations at the time of injection. Used as single dose prophylaxis, they should exert a negligible impact on resistance development. Their spectra encompass members of the aerobic and anaerobic (i.e. deep) skin flora, which may be introduced during the injection procedure, notably

staphylococci, hemolytic streptococci, and Cutibacterium (former Propionibacterium) spp. Dosages may need to be adjusted in patients with co-morbidities; hospital pharmacist advice should be sought, if necessary.

If these antibiotics are not available, please consult your national microbiological experts for recommendations.

PROCEDURE

The procedure must be conducted under aseptic conditions. Prior to injection, it is essential to swab at least 5 cm around the injection site using e.g., chlorhexidine with alcohol three times at one-minute intervals.

The injection of the hydrogel must be performed under local anaesthesia.

Injection of Arthrosamid® should be performed with the help of ultrasound in the lateral, proximal recess. If ultrasound is not available, the injection can be performed along the joint line lateral to the ligamentum patellae. The patient can be in a supine position or sitting with the knee bent between 80°- 90°.

Arthrosamid® is intended to be used with a sterile needle with Luer Lock fitting. The recommended needle is a standard 21G x 2 Inches (0.8 x 50 mm) with Luer Lock.

The procedure must be conducted using an aseptic no-touch technique. Insert the recommended needle into the knee joint cavity. If joint effusion is present, aspirate effusion and leave the needle in the knee.

Remove the protective Tip Cap from the Arthrosamid® syringe. Mount the Arthrosamid® syringe firmly onto the needle. Inject the hydrogel intra-articularly. Remove the syringe from the needle. Repeat steps for all syringes used for the procedure. In the event of gel leaking from the needle/syringe connection, secure the syringe more firmly. In the event of an obstruction of the needle, do not increase the pressure on the plunger. Instead, stop the injection and replace the needle. When the desired dosage has been administered, retract the needle from the injection site.

Recommended dosage of Arthrosamid® per treated knee: One single injection of 6 mL (6 x 1 mL).

Repeated injections of Arthrosamid® have not been evaluated in clinical studies.

A cumulative dose of more than 24 mL has not been evaluated in biological safety testing.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP is made available in the European database on medical devices (EUDAMED), where it is linked to the Basic UDI-DI.

URL to the EUDAMED public website: <https://ec.europa.eu/tools/eudamed>

Basic UDI-DI: 5704101Arthrosamid3K

The SSCP is available at <https://contura.com/orthopaedics/> until it is available in EUDAMED.

SYMBOLS USED ON LABELLING:

SYMBOL	EXPLANATION
	Manufacturer
	Sterile. Sterilized by moist heat
	Indication of Moist heat Sterile Barrier. Packaging level ensuring sterility
	Lot number
	Catalogue number
	Date of manufacture (YYYY-MM)
	Use before date printed on the label
	For single use only. Do not re-use
	Do not re-sterilize
	Keep away from sunlight
	Do not freeze
	Do not use if package is damaged and consult instructions for use
	Consult Instructions for Use
	Product is a medical device
	CE mark of conformity with applicable requirements set out in the Medical Device Regulations (EU) 2017/745
	Patient name or ID
	Implanting healthcare institution
	Website for patients
	UDI as AIDC format
	Date of implantation

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Tıbbi Cihaz Yönetmeliği (AB) 2017/745'te belirtilen geçerli gerekliliklere uygunluk CE işareti